

OptimAL Line 64 Bin Capacity Medication Cart, Standard Package #AL64MDBIN

OptimAL Line medication carts are lightweight due to the all-aluminum construction. These carts are available to accommodate all types of medication packaging and offer a complete line of accessories to meet your storage and security needs. Harloff carts are shipped fully assembled and ready to use.



Construction:

- Painted aluminum cabinet and drawers

Features:

- Internal narcotics box with Illinois Lock
- Ergonomic height
- BEST® brand lock on cabinet
- Painted pull-out side shelf with recessed work area
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with durable work surface
- Ergonomic top rail/push handle
- Easy clean protective pontoon bumpers
- 5" ball-bearing premium rubber casters, all full swivel, 1 brake and 1 directional
- Hard-baked powder coat finish
- Twelve year manufacturer's warranty

Drawer Configuration:

- 4 ea – 4" drawers
- 1 ea – 11" multipurpose storage drawer with locking narcotics box

Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>

Assembled Cart Dimensions:

- 42" H x 36" W x 24" D (includes top rail)
(106.6 cm x 91.4 cm x 61 cm)
- Writing surface: 38" (96.5 cm)
- Cart weight: 133 pounds (60.45 kg)
- Approx. shipping weight: 175 pounds (79.45 kg)