



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## Aluminum Universal Line Wide 4 Drawer Medication Cart for DISPILL® Packaging with CompX® Electronic Lock #AL8W10EKC4DSP

Harloff's Universal Line of medication carts were designed to accommodate all types of medication packaging. Its modern appearance incorporates industry-leading standard features and offers multiple locking options.

### Specifications

#### Construction:

- Aluminum, light-weight cabinet construction with 13" steel drawers

#### Features:

- Extended width single column medication cart with 100 to 150 DISPILL® card capacity
- Locking narcotics storage included
- Internal pullout shelf on right side
- CompX® brand electronic keypad lock on cabinet with key lock override
- Access rights may be assigned for up to 250 users per cart
- Drawers automatically re-lock after access (pre-set programmed time)
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumper
- Extended caster footprint for stability when drawers are open
- 5" ball-bearing premium casters, all full swivel- one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled and ready to use
- Pre-drilled hole locations for standard accessories
- Latex free
- Durable powder coat finish available in 17 colors
- Twelve year manufacturer's warranty

#### Internal Drawer Dimensions:

- 30" W x 17" D

#### Assembled Cart Dimensions:

- 48.5" H x 38.25" W x 25" D
- Writing Surface: 44"
- Cart weight: 145 pounds
- Approximate shipping weight: 180 pounds

#### Seventeen available color choices:

(may be viewed at [www.harloff.com](http://www.harloff.com))

Beige, Burgundy, Light Blue, Navy, Forest Green, Light Gray, Mauve, Red, Teal, Yellow, White, Hammer Black, Hammer Red, Hammertone Blue, Hammertone Gray, Hammertone Green, and Sand

#### Additional Option: 41052 – Lockview® Software

- USB dongle with RJ-11 cable.
- Allows user to interact with electronic lock through a laptop or PC.
- Software enables mass programming, view of audit trail (who accessed the cart and when), changes in auto re-lock timing, etc.

#### Drawer Configuration:

- (2) – 3.25" multipurpose storage drawer
- (1) – 13" drawer with locking narcotics box and rail system to hold DISPILL® packaging
- (1) – 13" drawer with rail system to hold DISPILL® packaging

