



MEDICATION CARTS **AL8W10K4DSP**

DISPILL® MEDICATION CART 132 Dispill® Card Capacity with a Key Lock

Features:

- Single column medication delivery cart for Dispill® cards and liquid storage with a maximum capacity of 132 Dispill® cards
- Includes one single width locking narcotics box made for Dispill® packaging that secures with a standard key lock
- Key locking cart with 4 keys per lock
- Durable scratch resistant heat sealed powder coat finish
- All aluminum construction
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Removeable plastic top
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Painted aluminum top rail / push handle
- Protective pontoon bumpers
- Extended caster footprint for stability when drawers are extended
- Body has pre-threaded accessory mounting holes
- Latex free cart is delivered fully assembled



Color shown is Cream. More than 20 additional colors are available.

Aluminum 132 Dispill® Card Capacity Medication Cart with Key Lock on Cabinet

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

Internal Drawer Dimensions:

- 30" W x 17" D

Assembled Cart Dimensions:

- Writing Surface: 44"
- Cart weight: 130 pounds (59 kg)
- 48.5" H x 38.25" W x 25" D (123.2 cm x 97.2 cm x 63.5 cm)
- Approx shipping weight: 170 pounds (77.1 kg)
- Approx shipping dims: 52" H x 46" W x 26" D (132.1 cm x 116.8 cm x 66 cm)

Drawer Configuration:

- 2 ea – 3.25" multipurpose storage drawer
- 1 ea – 13" drawer with locking narcotics box & rail system to hold Dispill® packs
- 1 ea – 13" drawer with rail system to hold Dispill® packs

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Engineered & Manufactured in the U.S.



View similar products online at
<https://www.harloff.com/product-category/specialty-packaging-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300