



## Slides Under Beds and Countertops

**ETC-3**

### ETC LINE TREATMENT CART with Three Drawers and a Lower Compartment

*Harloff's smallest cart, this economy treatment cart is shorter than standard counter heights and is designed to slide under equipment or beds.*

#### Features:

- Painted steel cabinet and drawers
- Low height make these carts ideal for bedside or under a counter
- Key lock with two keys
- 3" stem casters, two with brake
- Roller drawer slides
- Removable plastic top
- Latex free
- Heat sealed powder coat finish
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

Color shown is Red.  
More than 20 additional  
colors are available.



Engineered &  
Manufactured  
in the U.S.

### Three Drawer Treatment Cart with Lower Compartment



#### Drawer Configuration:

- (3) - 2.75-inch (6.98 cm)
- (1) - 12-inch (30.5 cm) Storage compartment

#### Assembled Cart Dimensions:

- Internal Drawer Dimensions: 20" W x 12.75" D (50.8 cm x 32.4 cm)
- 28.63" H x 23.5" W x 15" D (72.73 cm H x 59.4 W x 38.4 D)
- Cart weight: 48 lbs (21.8 kg)
- Approx. shipping weight: 76 pounds (34.5 kg)
- Approx. shipping dims: 29" H x 33" W x 23" D (73.7 cm H x 83.8 cm W x 58.4 cm D)

View similar products online at  
<https://www.harloff.com/m-series-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300

**HARLOFF**  
INNOVATING STORAGE FOR CARE