

ETC Line, Three Drawers with Lower Compartment, Treatment Cart #ETC-3



Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Low height make these carts ideal for bedside or under a counter
- Key lock with two keys
- 3" stem casters, two with brake
- Roller drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Twelve year warranty on manufactured parts
- Hard-baked powder coat finish

Drawer Configuration:

- (3) - 2.75-inch (6.98 cm)
- (1) - 12-inch (30.5 cm) Storage compartment

Internal Drawer Dimensions:

- 20" W x 12.75" D (50.8 cm x 32.4cm)

Assembled Cart Dimensions:

- 28.63" H x 23.5" W x 15" D
(72.73 cm H x 59.4 W x 38.4 D)
- Cart weight: 48 lbs (21.8 kg)
- Approx. shipping weight: 76 pounds (34.5 kg)
- Approx. shipping dims: 29" H x 33" W x 23" D
(73.7cm H x 83.8cm W x 58.4cm D)