



**M-series**

**M3DS1830K16**

## **M-SERIES TALL MINI CART**

Six Drawers with a Standard Key Lock

### **Features:**

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 30" of vertical drawer space
- Drawer Configuration: Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.4 cm)
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

### **Colors:**

Shown in Yellow, but available in more than twenty color options which you can view online at [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/).

All dimensions are available in the "MA-Dims-18wide-3DScaster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



**Engineered &  
Manufactured in  
the U.S.**



## **Tall Height, Mini Width, Six Drawer Infection Control Cart with Standard Key Lock**

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/product-category/m-series-medical-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**