

M-series

Mini Width Carts

HARLOFF

INNOVATING STORAGE FOR CARE

**Choose the Height and Width - Choose Your Drawer Configuration
Choose Your Lock - Your Color - Your Accessories**

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



M3DS1827K05



MDS1824EKC05



MDS1830B06



**Engineered &
Manufactured
in the U.S.**

Warranty:

- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

M-Series Features:

- Choice of lock options
- Full extension drawers slides for easy access
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Durable scratch resistant hard baked powder coat finish
- Industry best 5" quiet rolling Tente® casters
 - 1 steering caster & 2 braking casters
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable

Base MDS18 Cart Dimensions:

- Inside Drawer Dims: 12.5" W x 17" D (31.8 x 43.2 cm)
- Cart Weight – 85-105 lbs, depending on model
- Shipping Dims: 56" H x 29" W x 23" D (142.2 cm x 99 cm x 66 cm)

Toll Free: (800) 433-4064 • medicalsales@harloff.com • www.harloff.com • Ph: (719) 637-0300

The Choice is Yours!

X-Short
 34" H x 24.75" W x 22" D
 32.875" Height of Work Surface
 21" of Configurable
 Drawer Space



MDS1821E14

Short
 37.25" H x 24.75" W x 22" D
 36.125" Height of Work Surface
 24" of Configurable
 Drawer Space



M3DS1824K03

Medium
 40.5" H x 24.75" W x 22" D
 39.375" Height of Work Surface
 27" of Configurable
 Drawer Space




MDS1827WKC05

Tall
 43.75" H x 24.75" W x 22" D
 42.625" Height of Work Surface
 30" of Configurable Drawer Space



MDS1830EKC06



Choose Your Drawer Configuration:

1. Choose a cart size (X-Short, Short, Medium and Tall)
2. Select **any combination** of 3", 6", 9" or 12" drawers that fit within your cart's configurable drawer space.

Choose Your Lock:

- **B** - Breakaway Lock • **K** - Key-Lock • **E** - Basic Electronic Pushbutton Lock • **EKC** - Electronic Lock w/Keypad
 - **ELP** - Electronic Lock w/Keypad and Proximity Reader • **WKC** - Wireless Electronic Lock w/Keypad • **WELP** - Wireless Electronic Lock w/Keypad and Proximity Reader
- View locking options online: <https://www.harloff.com/locking-options/>

Choose Your Color(s):

Available in any of 23 colors. An anti-microbial finish is available in Lemon Drop Yellow. View color options online: <https://www.harloff.com/finish-options/>

Choose Your Accessories:

Full line of drawer dividers and exterior accessories are available for all cart sizes. View accessory options online: <https://www.harloff.com/accessories/>

Three M-Series Widths to Choose From:

