



M-series

MDS1830E06+MD18-ANS

M-SERIES TALL ANESTHESIA CART Mini Six Drawer Cart w/MD18-ANS Accessory Package

Tall Height, Mini Width, Six Drawer Anesthesia Cart with an Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD18-ANS Accessories Include:



3-Tier Back Rail (MD18-RBR3)



(3) 12" Side Rails (SIDERAIL12)



5-Compartment Tilt Bin w/Rail Clips (TILTBIN5RC)



Engineered & Manufactured in the U.S.

Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock with key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Cart is delivered assembled
- Fully customizable drawer configs with 30" of vertical drawer space
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Color shown with a Light Gray body and Hammertone Blue drawers. More than 20 additional colors are available.



View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Overall Cart Dimensions:

- Work Surface Height: 41.5" (105.4 cm)
- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Cart dims: 42.625" H x 23.875" W x 22" D (108.3 cm H x 60.6cm W x 55.9cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)