



M-series

MDS2421K05

M-SERIES MEDIUM WIDTH CART with Five Drawers and a Standard Key Lock

Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Ergonomic side push handles can be mounted on either side
- 21" of vertical drawer space
- Antimicrobial plastic top, removable
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - one steering caster and two braking casters
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Latex free



Engineered & Manufactured in the U.S.

Color shown with a Cream body and Navy drawers. More than 20 additional colors are available.



X-Short Height, Medium Width, Five Drawer Cart with Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Three 3" drawers (7.6 cm) and two 6" drawers (15.2 cm)

Assembled Cart Dimensions:

- 34" H x 30" W x 22" D (76.2 cm H x 76.2 cm W x 55.9 cm D)
- Cart net weight: 102 pounds (46.3 kg)
- Approx shipping weight: 141 pounds (64 kg)
- Approx shipping dims: 50" H x 33" W x 26" D (127 cm x 99 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300