



# M-series

# MDS3030K06+MD30-ANS

## M-SERIES TALL ANESTHESIA CART with MD30-ANS Anesthesia Package

### Tall, Standard Width, Six Drawer Anesthesia Cart with Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered &  
Manufactured  
in the U.S.

### MD30-ANS Accessories Include:



3-Tier Back Rail  
(MD30-RBR3)



(3) 12" Side Rails  
(SIDERAIL12)



5-Compartment  
Tilt Bin w/Rail Clips  
(TILTBIN5RC)



### Features:

- Sturdy, uni-body steel construction
- Choice of lock options: [www.harloff.com/locking-options/](http://www.harloff.com/locking-options/)
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Ergonomic side push handles can be mounted on either side
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Antimicrobial plastic top, removable
- Fully customizable drawer configurations
- 30" of vertical drawer space
- Cart is delivered assembled
- Latex free

Color shown is Navy.  
More than 20 additional  
colors are available.

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)
- 24.5" W x 17" D (62.2 x 43.2 cm)

#### Assembled Cart Dimensions:

- 66.75" H x 37.5" W x 22" D (169.5 cm H x 95.3 cm W x 55.9 cm D)
- Cart weight: 198 lbs (89.8 kg)
- Approximate shipping weight: 237 lbs (107.5 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

# HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300  
[medicalsales@harloff.com](mailto:medicalsales@harloff.com) • [www.harloff.com](http://www.harloff.com)