

MALIGNANT Hyperthermia

MH5100K

MALIGNANT HYPERTHERMIA CART

with 1.8 Cubic Feet Medical Grade Refrigerator

Features:

- 20-gauge powder-coated steel construction w/built-in Follett® REF2 refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- Highest quality in the industry Follett® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawer for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- · Key locks for drawer and refrigerator
- · Cart ships ready to use with wheels installed
- · Hospital grade power strip and cord included
- · Antimicrobial plastic top, removable


Colors:

Warranty:

Shown in Sand, but available in over twenty color options. View them online: www.harloff.com/finish-options/

Malignant Hyperthermia Cart, One Drawer with Key Lock

Accessories Included:





Waste container &

mounting bracket



document holders



Engineered & Manufactured in the U.S. Two large chart/

- Approximate shipping weight: 215 lbs (97.7 kg) Large wire basket
- Approximate shipping dims: 52" H x 44" W x 30" D (132.1cm H x 111.8cm W x 76.2cm D)

12-Year Manufacturer's Warranty on shipments

Drawer Configuration:

Overall Cart Dimensions: Internal Drawer Dimensions:

(1) 3" (7.6 cm) and (1) Compartment with pre-installed refrigerator

Cart dims: 45.5" H x 37.5" W x 26" D

(115.6 cm H x 95.3cm W x 66cm D) Cart net weight: 175 lbs (79.5 kg)

OPTIONAL:

Side storage cabinets with key locking door and two adjustable shelves for additional supply storage (each side cabinet replaces existing side mounted accessories).



View similar products online at

https://www.harloff.com/malignant-hyperthermia-cart/

Toll Free: (800) 433-4064 • medicalsales@harloff.com www.harloff.com • Ph: (719) 637-0300

