



# MALIGNANT HYPERTHERMIA



## MH5200K-AC

### MALIGNANT HYPERTHERMIA CART

with 2.4 Cubic Feet Medical Grade Refrigerator

#### Features:

- 20-gauge powder-coated steel construction w/built-in Accucold FF28LWHMED2 refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- High quality Accucold refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Key locks for drawers and refrigerator
- Cart ships ready to use with wheels installed
- Antimicrobial plastic top, removable

#### Colors:

Shown in Hammertone Blue,

but available in over twenty color options.

View them online: [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/)



### Malignant Hyperthermia Cart, Two Drawers with Key Lock

#### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

#### Overall Cart Dimensions:

- Internal Drawer Dimensions:  
23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 48.5" H x 37.5" W x 26" D  
(123.2 cm H x 95.3cm W x 66cm D)
- Cart net weight: 160 lbs (72.5 kg)
- Approximate shipping weight: 199 lbs (90.3 kg)
- Approximate shipping dims: 62" H x 44" W x 30" D  
(157.5cm H x 111.8cm W x 76.2cm D)

#### Drawer Configuration:

- (2) 3" (7.6 cm) and (1) Compartment with pre-installed refrigerator

#### Accessories Included:



Large wire basket



Waste container & mounting bracket



Two large chart/document holders



Engineered &  
Manufactured  
in the U.S.

#### OPTIONAL:

Side storage cabinets with key locking door and two adjustable shelves for additional supply storage (each side cabinet replaces existing side mounted accessories).



View similar products online at

<https://www.harloff.com/product-category/malignant-hyperthermia-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)

[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300