



# MALIGNANT HYPERTHERMIA

## MH5216B

### MALIGNANT HYPERTHERMIA CART

with 1.0 Cubic Feet Medical Grade Refrigerator

#### Features:

- 20-gauge powder-coated steel construction w/built-in Follett® REF1 refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- Highest quality in the industry Follett® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Breakaway locks on drawer and refrigerator (refrigerator can be keylocked)
- Cart ships ready to use with wheels installed
- Removable plastic top

#### Colors:

Shown in Navy, but available in over twenty color options.

View them online: [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/)



### Malignant Hyperthermia Cart, Two Drawer with Breakaway Lock

#### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

#### Overall Cart Dimensions:

- Internal Drawer Dimensions:  
23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 45.5" H x 37.5" W x 26" D  
(115.6 cm H x 95.3 cm W x 66 cm D)
- Cart net weight: 175 lbs (79.5 kg)
- Approximate shipping weight: 215 lbs (97.7 kg)
- Approximate shipping dims: 52" H x 44" W x 30" D  
(132.1 cm H x 111.8 cm W x 76.2 cm D)

#### Drawer Configuration:

- (2) 6" (15.2 cm) and (1) Compartment with pre-installed refrigerator

#### Accessories Included:



Large wire basket



Waste container & mounting bracket



Two large chart/document holders



Engineered &  
Manufactured  
in the U.S.

#### OPTIONAL:

Side storage cabinets with key locking door and two adjustable shelves for additional supply storage (each side cabinet replaces existing side mounted accessories).

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/malignant-hyperthermia-cart/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)

[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300