



# MALIGNANT HYPERTHERMIA

## MH5216K

### MALIGNANT HYPERTHERMIA CART with 1.0 Cubic Feet Medical Grade Refrigerator

#### Features:

- 20-gauge powder-coated steel construction w/built-in Follett® REF1 refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- Highest quality in the industry Follett® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Key locks for drawer and refrigerator
- Cart ships ready to use with wheels installed
- Hospital grade power strip and cord included
- Antimicrobial plastic top, removable

#### Colors:

Shown in Navy, but available in over twenty color options. View them online: [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/)



### Malignant Hyperthermia Cart, Two Drawer with Key Lock

#### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

#### Overall Cart Dimensions:

- Internal Drawer Dimensions: 23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 45.5" H x 37.5" W x 26" D (115.6 cm H x 95.3cm W x 66cm D)
- Cart net weight: 175 lbs (79.5 kg)
- Approximate shipping weight: 215 lbs (97.7 kg)
- Approximate shipping dims: 52" H x 44" W x 30" D (132.1cm H x 111.8cm W x 76.2cm D)

#### Drawer Configuration:

- (2) 6" (15.2 cm) and (1) Compartment with pre-installed refrigerator

#### Accessories Included:



Large wire basket



Waste container & mounting bracket



Two large chart/document holders



Engineered & Manufactured in the U.S.

#### OPTIONAL:

Side storage cabinets with key locking door and two adjustable shelves for additional supply storage (each side cabinet replaces existing side mounted accessories).



View similar products online at <https://www.harloff.com/malignant-hyperthermia-cart/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300