



MALIGNANT HYPERTHERMIA

MH5300B

MALIGNANT HYPERTHERMIA CART

with 1.8 Cubic Feet Medical Grade Refrigerator

Features:

- 20-gauge powder-coated steel construction w/built-in Follett® REF2 refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- Highest quality in the industry Follett® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Breakaway locks on drawer and refrigerator (refrigerator can be keylocked)
- Cart ships ready to use with wheels installed
- Antimicrobial plastic top, removable

Colors:

Shown in Beige, but available in over twenty color options.

View them online: www.harloff.com/finish-options/



Malignant Hyperthermia Cart, Three Drawers with Breakaway Lock

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Overall Cart Dimensions:

- Internal Drawer Dimensions:
23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 51.5" H x 37.5" W x 26" D
(130.8 cm H x 95.3cm W x 66cm D)
- Cart net weight: 188 lbs (85.3 kg)
- Approximate shipping weight: 299 lbs (135.6 kg)
- Approximate shipping dims: 62" H x 44" W x 30" D
(157.5cm H x 111.8cm W x 76.2cm D)

Drawer Configuration:

- (3) 3" (7.6 cm) and (1) Compartment with pre-installed refrigerator

Accessories Included:



Large wire basket



Waste container & mounting bracket



Two large chart/document holders



Engineered &
Manufactured
in the U.S.

OPTIONAL:

Side storage cabinets with key locking door and two adjustable shelves for additional supply storage (each side cabinet replaces existing side mounted accessories).

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/malignant-hyperthermia-cart/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300