

# MALIGNANT HYPERTHERMIA

# MH5300B

# **MALIGNANT HYPERTHERMIA CART**

with 1.8 Cubic Feet Medical Grade Refrigerator

#### Features:

- 20-gauge powder-coated steel construction w/built-in Follett® REF2 refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- Highest quality in the industry Follett® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Breakaway locks on drawer and refrigerator (refrigerator can be keylocked)
- · Cart ships ready to use with wheels installed
- · Hospital grade power strip and cord included
- · Antimicrobial plastic top, removable

#### Colors:

Shown in Hammertone Blue, but available in over twenty color options. View them online: <a href="www.harloff.com/finish-options/">www.harloff.com/finish-options/</a>

## Malignant Hyperthermia Cart, Three Drawers with Breakaway Lock



#### Warranty:

 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other Internationalshipments.
 1-Year Limited Warranty on any bought-out parts

#### **Overall Cart Dimensions:**


- Internal Drawer Dimensions:
  23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 51.5" H x 37.5" W x 26" D (130.8 cm H x 95.3cm W x 66cm D)
- Cart net weight: 188 lbs (85.3 kg)
- Approximate shipping weight: 227 lbs (103.2 kg)
- Approximate shipping dims: 62" H x 44" W x 30" D (157.5cm H x 111.8cm W x 76.2cm D)

#### **Drawer Configuration:**

• (3) 3" (7.6 cm) and (1) Compartment with pre-installed refrigerator

### **Accessories Included:**







Waste container & mounting bracket



Two large chart/ document holders



Engineered & Manufactured in the U.S.

#### **OPTIONAL:**

Side storage cabinets with key locking door and two adjustable shelves for additional supply storage (each side cabinet replaces existing side mounted accessories).

#### View similar products online at

https://www.harloff.com/malignant-hyperthermia-cart/

Toll Free: (800) 433-4064 • medicalsales@harloff.com www.harloff.com • Ph: (719) 637-0300

