



MPA1830EKC06

ALUMINUM

A-SERIES MINI WIDTH CART w/Six Drawers and an Electronic Lock w/Keypad

Features:

- Extremely light weight cart with critical grade aluminum construction, ideal for carts that are moved consistently throughout the day
- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Removable plastic top
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems. View online at <https://www.harloff.com/locking-options/>
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.9 cm)
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Color shown in
Hammerstone
Blue. More
than 20
additional
colors are
available.



**Engineered &
Manufactured
in the U.S.**

Tall Height, Mini Width, Six Drawer Cart with an Electronic Lock with Keypad

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.



View similar products online at
<https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300