



MPA1830ELP06


ALUMINUM

A-SERIES MINI WIDTH CART
with Six Drawers, Electronic Keypad Lock & Prox

Features:

- Extremely light weight cart with critical grade aluminum construction, ideal for carts that are moved consistently throughout the day
- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Antimicrobial plastic top, easy to remove/replace
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling Tente® casters - one steering caster and two braking casters
- Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems. View online at <https://www.harloff.com/locking-options/>
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space

Engineered & Manufactured in the U.S.



Color shown in Navy with Light Blue drawers. More than 20 additional colors are available.

Tall Height, Mini Width, Six Drawer Cart and Electronic Lock with Keypad and Prox Reader

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.9 cm)
- 12.5" W x 17" D (31.8 cm x 43.2 cm)

Assembled Cart Dimensions:

- 43.75" H x 23.875" W x 22" D (111.1 cm H x 60.6 cm W x 55.9 cm D)
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm x 61 cm x 61 cm)

Warranty:

- 1-Year Limited Warranty on any bought-out parts.
- 12-Year Manufacturer's Warranty on shipments to US and Canada.



View similar products online at <https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300