



# MPA3027K06+MD30-ANS

## ALUMINUM

### A-SERIES STANDARD WIDTH CART w/Six Drawers, Key Lock & MD30-ANS Anesthesia Pkg

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.

#### MD30-ANS Accessories Include:



3-Tier Back Rail  
(MD30-RBR3)



(3) 12" Side Rails  
(SIDERAIL12)



5-Compartment  
Tilt Bin w/Rail Clips  
(TILTBIN5RC)



Engineered &  
Manufactured  
in the U.S.



- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Antimicrobial plastic top, removable and easily replaceable
- Fully customizable drawer configurations at no extra charge
- Ergonomic side push handles, mounted on the left or right
- Industry best 5" quiet rolling Tente® casters
  - 1 steeringcaster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space

Color shown is Navy.  
More than 20 additional  
colors are available.

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)
- 24.5" W x 17" D (62.2 x 43.2 cm)

#### Assembled Cart Dimensions:

- 66.75" H x 37.5" W x 22" D (169.5 cm H x 95.3 cm W x 55.9 cm D)
- Cart weight: 198 lbs (89.8 kg)
- Approximate shipping weight: 237 lbs (107.5 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300  
medicalsales@harloff.com • www.harloff.com