



MRI-SAFE



MR6B

MR-CONDITIONAL EMERGENCY CART with Six Drawers and a Breakaway Lock

Features:

- Tested and certified by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; two swivel – two with brake, one directional
- Breakaway locking system (with optional multiple lock bars)
- Pack of 100 breakaway seals
- Protective corner bumpers
- Pull-out shelf
- Removable plastic top with integrated push handles
- Latex free cart is delivered fully assembled with pre-threaded accessory mounting holes
- Durable, heat sealed powder coat finish



Engineered &
Manufactured
in the U.S.



Color shown
is White with
Red drawers.
More than
20 additional
colors are
available.

Drawer Configuration:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)
- One 4" drawer (10.2 cm), two 3.25" drawers (8.3 cm), two 6.5" drawers (16.5 cm) and one 9.75" drawer (24.8 cm)

Carts Specially Made for Magnetic Resonance Imaging Suites

Harloff offers a line of carts that meet American Society for Testing and Materials (ASTM) International Designation of MR-Conditional for use in a 3 Tesla or less environment of the Magnetic Resonance Imaging Suite (testing date 1-13-2010). Whether you need it for Anesthesia or Emergency, Harloff has the right MR-Conditional cart for your needs.

Assembled Cart Dimensions:

- Writing surface: 43.5" (110.5cm)
- 46" H x 31.8" W x 23.75" D (116.8 cm H x 80.8 cm W x 60.3 cm D)
- Cart net weight: 108 lbs (49.5 kg)
- Approximate shipping weight: 146 lbs (66.4 kg)
- Approximate shipping dims: 55" H x 37" W x 26" D (139.7cm H x 94cm W x 66cm D)

Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/mri-safe-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300