

MR-Conditional Line/ Model: MR6B

Aluminum Six Drawer Anesthesia Cart, Breakaway Lock, Standard Package



Over 20 Drawer Colors Available:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Burgundy, Red, Yellow, Beige, Sand, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Cream, Taupe, White, Brushed Silver, Forest Green, HammerTone Black, HammerTone Red and Pink.

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Non-magnetic or "weekly magnetic" materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

Features:

- Designed for compatibility in MRI suites
- Tested by independent testing authority to Tesla 3 level (copies of report available)
- Meets American Society for testing and materials (ASTM) international, designation: F2503-5 de inition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; two swivel – two with brake, one directional
- Painted aluminum break-away lock bar with 100 break-away seals
- Corner bumpers, protecting the cart
- Painted aluminum pull-out shelf
- Antimicrobial plastic top, removable
- MR-Conditional accessories available
- Latex free
- Durable, hard-baked powder coat inish

Drawer Configuration:

- (1) - 4" (10.2 cm)
- (2) - 3.25" (8.3 cm)
- (2) - 6.5" (16.5 cm)
- (1) - 9.75" (24.8 cm)

Internal Drawer Configuration:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)

Overall Cart Dimensions:

- 45.9" H x 30.5" W x 23.9" D (116.2 cm H x 77.5 cm W x 66.7 cm D)
- Writing surface: 43.5" (110.5 cm)
- Cart net weight: 109 pounds (49.5 kg)
- Approximate shipping weight: 148 pounds (67.3 kg)
- Approximate shipping dims: 55" H x 37" W x 26" D (139.7 cm H x 94 cm W x 66 cm D)