



MRI-SAFE

MR7K

MR-CONDITIONAL ANESTHESIA CART with Seven Drawers and a Key Lock

Features:

- Tested and certified by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Protective corner bumpers
- Pull-out shelf
- Low-ferrous 5" sealed ball-bearing casters; two swivel – two with brake, one directional
- BEST® brand lock on cabinet, with 2 keys
- Removable plastic top with integrated push handles
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- MR-Conditional accessories available
- Latex free
- Durable, heat sealed powder coat finish



**Engineered &
Manufactured
in the U.S.**



**Color shown
is White with
Hammertone
Green drawers.
More than 20
additional
colors are
available.**

Drawer Configuration:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)
- One 4" drawer (10.2 cm), four 3.25" drawers (8.3 cm), one 6.5" drawer (16.5 cm) and one 9.75" drawer (24.8 cm)

Carts Specially Made for Magnetic Resonance Imaging Suites

Harloff offers a line of carts that meet American Society for Testing and Materials (ASTM) International Designation of MR-Conditional for use in a 3 Tesla or less environment of the Magnetic Resonance Imaging Suite (testing date 1-13-2010). Whether you need it for Anesthesia or Emergency, Harloff has the right MR-Conditional cart for your needs.

Assembled Cart Dimensions:

- Writing surface: 43.5" (110.5cm)
- 46" H x 31.8" W x 23.75" D (116.8 cm H x 80.8 cm W x 60.3 cm D)
- Cart net weight: 137 lbs (62.3 kg)
- Approximate shipping weight: 157 lbs (71.4 kg)
- Approximate shipping dims: 55" H x 37" W x 26" D (139.7cm H x 94cm W x 66cm D)

Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/product-category/mri-safe-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300