

## MR-Conditional Line/ Model: MRN6B

Narrow Aluminum Six Drawer Emergency Cart, Breakaway Lock, Standard Package



### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

### Over Twenty Available Drawer Colors:

**Antimicrobial Finish Available: Lemon Drop Yellow Only**  
Burgundy, Red, Yellow, Beige, Sand, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Cream, Taupe, White, Brushed Silver, Forest Green, HammerTone Black, HammerTone Red and Pink.

### Construction:

- Non-magnetic or "weekly magnetic" materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

### Features:

- Designed for compatibility in MRI suites
- Tested by independent testing authority to 3 Tesla level (copies of report available)
- Meets American Society for Testing and Materials (ASTM) International, designation: F2503-05 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel – one with brake, one directional
- BEST® brand lock on cabinet, with 2 keys
- Painted aluminum break-away lock bar with 100 break-away seals
- Corner bumpers, protecting the cart
- Antimicrobial plastic top, removable
- MR-Conditional accessories available
- Latex free
- Durable, hard-baked powder coat finish

### Drawer Configuration:

- (1) - 4" (10.2 cm)
- (2) - 3.25" (8.3 cm)
- (2) - 6.5" (16.5 cm)
- (1) - 9.75" (24.8 cm)

### Internal Drawer Configuration:

- 13.75" W x 17" D (34.9 cm x 88.7 cm)

### Overall Cart Dimensions:

- 46" H x 20" W x 23.5" D (116.8 cm H x 50.8 cm W x 59.7 cm D)
- Writing surface: 44" (111.7 cm)
- Cart net weight: 85 lbs (38.6 kg)
- Approximate shipping weight: 123 lbs (55.9 kg)
- Approximate shipping dims: 57" H x 33" W x 26" D (144.8 cm H x 83.8 cm W x 66 cm D)