

MR-Conditional Line/ Model: MRN6K-MAN

Narrow Aluminum Six Drawer Anesthesia Cart, Key Lock, Specialty Package



Construction:

- Non-magnetic or “weekly magnetic” materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

Warranty:

- 12-Year Manufacturer’s Warranty on shipments to US and Canada, 5-Year Manufacturer’s Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Over 20 Drawer Colors Available:

Antimicrobial Finish Available: Lemon Drop Yellow Only

Burgundy, Red, Yellow, Beige, Sand, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Cream, Taupe, White, Brushed Silver, Forest Green, HammerTone Black, HammerTone Red and Pink.

Features:

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of “MR-Conditional”
- Narrow design for use in MRI suite with limited floor space
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5” sealed ball-bearing casters; full swivel – one with brake, one directional
- BEST® brand lock on cabinet, with 2 keys
- Corner bumpers, protecting the cart
- Replaceable plastic top with integrated push handles
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- MR-Conditional accessories available
- Latex free
- Hard-baked powder coat finish

Anesthesia Specialty Package Includes MRI-Compatible Accessories:

- Aluminum 2-tier raised back rail system (MR-RBR2)
- Aluminum side rails (4 ea.) (MR-SR)
- Plastic full drawer tray organizers (2 ea.) (68530-P1)

Drawer Configuration:

- (1) - 4” (10.2 cm)
- (2) - 3.25” (8.3 cm)
- (2) - 6.5” (16.5 cm)
- (1) - 9.75” (24.8 cm)

Internal Drawer Configuration:

- 13.75” W x 17” D (34.9 cm x 88.7 cm)

Overall Cart Dimensions:

- 68.5” H x 20” W x 23.5” D (174 cm H x 50.8 cm W x 59.7 cm D)
- Writing surface: 44” (111.7 cm)
- Cart net weight: 108 lbs (50.5 kg)
- Approximate shipping weight: 144 lbs (65.5 kg)
- Approximate shipping dims: 57” H x 33” W x 26” D (144.8 cm H x 83.8 cm W x 66 cm D)