

MR-Conditional Line/ Model: MRN6K

Narrow Aluminum Six Drawer Anesthesia Cart, Key Lock, Standard Package



Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Over Twenty Drawer Colors Available:

Antimicrobial Finish Available: Lemon Drop Yellow Only

Burgundy, Red, Yellow, Beige, Sand, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Cream, Taupe, White, Brushed Silver, Forest Green, HammerTone Black, HammerTone Red and Pink.

Construction:

- Non-magnetic or "weakly magnetic" materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

Features:

- Designed for compatibility in MRI suites
- Tested by independent testing authority to Tesla 3 level (copies of report available)
- Meets American Society for testing and materials (ASTM) international, designation: F2503-5 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel – one with brake, one directional
- BEST® brand lock on cabinet, with 2 keys
- Corner bumpers, protecting the cart
- Antimicrobial plastic top, removable
- MR-Conditional accessories available
- Latex free
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (1) - 4" (10.2 cm)
- (2) - 3.25" (8.3 cm)
- (2) - 6.5" (16.5 cm)
- (1) - 9.75" (24.8 cm)

Internal Drawer Configuration:

- 13.75" W x 17" D (34.9 cm x 88.7 cm)

Overall Cart Dimensions:

- 46" H x 20" W x 23.5" D (116.8 cm H x 50.8 cm W x 59.7 cm D)
- Writing surface: 44" (111.7 cm)
- Cart net weight: 90 lbs (40.9 kg)
- Approximate shipping weight: 123 lbs (55.79 kg)
- Approximate shipping dims: 57" H x 33" W x 26" D (144.8 cm H x 83.8 cm W x 66 cm D)