



MRI-SAFE

MRN7B-EMG

MR-CONDITIONAL EMERGENCY CART 7 Drawers, Breakaway Lock, Specialty Package


Specialty Package Includes:



Oxygen tank holder (MR-ALO2)



Cardiac board and brackets (MR-CARDBRD)



(2) Plastic Trays w/ Adjustable Dividers (68620-P1)



(4) Utility hooks (MR-UHOOKS4)



IV pole (MR-IV2)

Features:

- Tested and certified by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; two swivel – two with brake, one directional
- Breakaway locking system (with optional multiple lock bars)
- Pack of 100 breakaway seals
- Protective corner bumpers
- Removable plastic top with integrated push handles
- Latex free cart is delivered fully assembled with pre-threaded accessory mounting holes
- Durable, heat sealed powder coat finish



Engineered & Manufactured in the U.S.

Color shown with a White body and Red drawers. More than 20 additional colors are available.



Drawer Configuration:

- 13.75" W x 17" D (34.9 cm x 88.7 cm)
- One 4" drawer (10.2 cm), two 3.25" drawers (8.3 cm), two 6.5" drawers (16.5 cm) and one 9.75" drawer (24.8 cm)

Assembled Cart Dimensions:

- Writing surface: 44" (111.7 cm)
- 46" H x 21.5" W x 23.75" D (116.8 cm H x 54.6 cm W x 60.3 cm D)
- Cart net weight: 116 lbs (52.7 kg)
- Approximate shipping weight: 152 lbs (69.1 kg)
- Approximate shipping dims: 57" H x 33" W x 26" D (144.8 cm H x 83.8 cm W x 66 cm D)

Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Specially Made for Magnetic Resonance Imaging Suites



View similar products online at <https://www.harloff.com/product-category/mri-safe-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300