



- 840 Box Capacity
- Mobile Body
- Glass Doors
- Standard Key Lock

**MEDSTOR<sup>MAX</sup>**

**MS-IOL840**

## INTRAOCULAR LENS CART

### Features:

- Unmatched quality of construction — built to last using professional grade design and manufacturing for medical-grade applications
- Key lock with two keys for lockable, secure storage of supplies
- Features a true 24" interior cabinet depth for storage of up to 840 IOL boxes
- Four pull-out drawer shelves with adjustable dividers — ideal for medical supplies used in surgery centers and operating rooms
- Double doors can swing up to 270 degrees for easy access
- Glass windows allow for visibility of supplies and inventory
- Ball-bearing full extension drawer slides
- Industry best 5" quiet rolling casters for superior mobility - one steering and one braking caster
- Cabinet is delivered fully assembled

**Color shown in Hammertone Blue. More than 20 additional colors are available.**

#### Assembled Cart Dimensions:

- 40.75" H x 40.1" W x 29.25" D  
(103.5 cm H x 101.9 cm W x 74.3 cm D)
- Cart net weight: 256 pounds (116.1 kg)
- Approximate shipping weight: 400 lbs (181.4 kg)
- Approximate Shipping Dims: 62" H x 49" W x 39" D (157.5 cm x 124.5 cm x 99.1 cm)

**Warranty:** 12-Year Limited Manufacturer's Warranty,  
1-Year Limited Warranty on any bought-out parts



### Mobile IntraOcular Lens Storage Cart, Half Height Cabinet, Four Drawers, Glass Doors, Standard Key Lock

Attractive, secure storage for IntraOcular lenses. Provides mobility and is compatible with multiple brands and sizes of IOL boxes including Alcon, Bausch and Lomb, Rayner, and Abbot Medical Optics.



**Engineered & Manufactured in the U.S.**

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/product-category/intraocular-lens-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**