



**Integrated Audio  
and Visual alarm  
that alerts when open**

**NC16C12-DT2-AVD**

**MEDIUM NARCOTICS CABINET**  
with Locking Outer Door & Locking Inner Door

**Construction:**

- Heavy duty 20 gauge steel


**Features:**

- Cabinet has integrated Audio and Visual alarm to alert personnel when open
- Battery-powered using standard 9 volt; battery in quick-change holder
- Does not need to be hard-wired to electrical power - can be mounted anywhere
- Dual inner door and outer door construction
- Two high security pick-resistant tubular locks
- Four keys per lock included
- Unique lock sequences– all cabinets are keyed differently
- Full length stainless steel pinned door hinges prevent pin removal
- One adjustable and/or removable shelf with the choice of five locations
- Predrilled holes in the back for mounting on a wall, cabinet or other solid surface (wall mounting hardware included)
- Narcotics cabinet is delivered ready to use
- Durable heat sealed powder coat finish



**Engineered & Manufactured in the U.S.**

*This narcotics cabinet fully meets DEA guidelines for the secure storage of Schedule I and II controlled substances provided that the cabinet is bolted to a wall, floor or other solid, non-mobile surface and provided that other guidelines are followed.*



**Shown in gloss beige (default color). More than 20 additional colors are available.**

**Keys with black tips must be used in locks with black ring. Putting the wrong key in a lock will disable the lock and require the lock to be replaced.**



**One lock and key set includes a black identifier**

Harloff offers many different narcotic cabinets to help you keep both your controlled substances and your business secure including wood look cabinets, slim profile boxes and cabinets with audio/visual alarms.

**Assembled Cart Dimensions:**

- 16" H x 12" W x 8" D (41 cm H x 30.5 cm W x 20.4 cm D)
- Cabinet weight: 15 lbs (6.8 kg)
- Approximate shipping weight: 17 lbs (7.7 kg)
- Approximate shipping dims: 18" H x 16" W x 12" D (45.7 cm H x 40.6 cm x 30.5 cm D)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts



**View similar products online at**  
<https://www.harloff.com/product-category/narcotics-medicine-cabinets/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**