



**For punch cards,  
pill bottles and  
liquid storage.**

# NC30D30-BK2



**LARGE NARCOTICS CABINET**  
with a Double Door and Two Key Locks

## Features:

- Heavy duty 20 gauge steel construction
- Double locking cabinet, double door
- Each shelf holds 4 rows of 6" x 9" punch cards with a total capacity of 180 punch cards
- Open area at bottom for liquid storage and misc supplies
- Two keys per lock; master keys are available
- Full length stainless steel pinned door hinge prevents pin removal
- One fixed shelf and one adjustable/removable shelf with the choice of 5 shelf locations
- Predrilled holes in the back for mounting on a wall, cabinet or other solid surface (wall mounting hardware included)
- Narcotics cabinet is delivered ready to use
- Durable heat sealed powder coat finish

**Engineered &  
Manufactured  
in the U.S.**



**Contents  
not  
included.**



**Shown in gloss beige  
(default color).  
More than 20 additional  
colors are available.**

*This cabinet is compliant with DEA guidelines for the secure storage of Schedule III, IV, and V controlled substances when a separate outer door is also used, such as inside a drawer or cabinet within a locked office or in a locked refrigerator/freezer within a locked office. Because this narcotics cabinet does not have both an inner and outer door, it is not recommended for the storage of Schedule I and II controlled substances unless additional criteria is met. Please consult current DEA guidelines for details.*

### Assembled Cart Dimensions:

- 30" H x 30" W x 10" D (76.2 cm H x 76.2 cm W x 25.4 cm D)
- Net weight: 54 pounds (24.5 kg)
- Approx. shipping weight: 62 pounds (28.1 kg)
- Approximate shipping dims: 34" H x 34" W x 14" D (86.4 cm x 86.4 cm x 35.6 cm)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

## Large Narcotics Cabinet with a Double Door and Two Key Locks

Harloff offers many different narcotic cabinets to help you keep both your controlled substances and your business secure including wood look cabinets, slim profile boxes and cabinets with audio/visual alarms.

**View similar products online at**

<https://www.harloff.com/product-category/narcotics-medicine-cabinets/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**

