

FOR BRONCHOSCOPES

ENT scopes & other small scopes



ureDry
SCOPE CABINETS

**FOUR BRONCHOSCOPE OR EIGHT
ENT SCOPE CAPACITY WALL MOUNT
STORAGE CABINET, KEY LOCK**

Scopes not included

SCW2448DR

*This new wall mounted
bronchoscope cabinet
provides the same safe,
vertical drying capabilities
as that of scopes in other
parts of your facility*

- Customizable scope holders to fit most common scope types
- Holds up to 4 Bronchoscopes or up to 8 ENT Scopes
- Glass paneled doors for easy viewing of scope inventory
- Air flow vents on the side of the cabinet
- Industry best 12 year limited manufacturer warranty
- Cabinet dimensions:
48" H x 24" W x 12.75" D



**Engineered &
Manufactured in the U.S.**

*Comes in Stainless Steel
or heat sealed powder
coat steel, available in
over 20 color options*

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/ent-and-bronchoscope-cabinets/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**



10 Reasons to Choose **HARLOFF** for Your Scope Cabinets

A Smart Choice to Meet & Exceed Current SGNA, AORN and AAMI Guidelines



1

Industry Best Positive Pressure HEPA Filtered Drying System. Designed to add Dri-Scope Aid® now or later with easy “Plug & Play” installation.

Proudly Made in the USA



Harloff products are proudly **Designed, Engineered and Manufactured in the United States** at our Colorado Springs, CO facility.

6

Drying or Non-Drying Options

For basic storage and air drying of reprocessed scopes, this model is the perfect solution. **Need HEPA filtered positive pressure scope drying too?** Upgrade to the DRDP version for the quickest removal of biofilm-producing moisture using filtered air.

2

Custom Requests - No Problem!

Our in-house engineering team gives us an unparalleled ability to do custom and unique projects. *Let us design the perfect scope storage system for your facility!*

7

Optional Wall or Cart Mounting

Our cabinets fit your floor plan! Use our wall-mounted cabinet for permanent installations or choose a mobile cart-mounted cabinet to move scopes from your reprocessing area directly to your treatment room.

3

Keyed or Keyless Locks

Options include latch locks, key locks or standard electronic push button locks. View

your locking options online at www.harloff.com/locking-options/

8



Versatile Scope Holders

Holds major brands including **Optim, Fuji, Olympus, Pentax, Karl Storz and others.**

4

Scope Cabinet Finish Options

Choose Stainless Steel, Heat Sealed Powder Coat in High Gloss or Textured Low Gloss, or upgrade to an **antimicrobial finish.**

9

12 Year Warranty

Our Manufacturer's Warranty is Unmatched in the industry. We build our cabinets to last and we stand behind them.

5



Superior Color Choices

Over 20 Colors to choose from. Cabinets have a heat sealed powder coat finish inside and out. View color options online at www.harloff.com/finish-options/

10