

Standard Line/ Model: SL270MOT

Medication Cart, key Lock, Standard Package



Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Single column medication delivery cart
- Capacity of 270 Medicine On Time® cards
- Key lock system with Best® brand lock on cabinet
- Illinois lock on narcotics drawer located in top 9" drawer, designed to hold MOT® cards
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- 5" ball-bearing premium casters, all full swivel, 2 with toe locking brake & 1 directional
- Pre-threaded accessory mounting holes
- Corner bumpers, protecting the cart

Internal Drawer Configuration:

- 17" D x 29.5" W

Drawer Configuration:

- (1) – 4" multipurpose storage drawer
- (1) – 9.75" punch card drawer with dividers and locking narcotics box
- (2) – 9.75" punch card drawers with dividers

Overall Cart Dimensions:

- 45.75" H x 36" W x 23.5" D (116.2 cm H x 91.4 cm W x 59.7 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 165 lbs (75 kg)
- Approx. shipping weight: 220 lbs (100 kg)
- Approx. shipping dims: 57" H x 36" W x 26" D (144.8 cm H x 99.1 cm W x 66 cm D)