



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Standard Line Multi-Dose Bin Medication Cart, Standard Package #SL72MD

### Specifications

#### Construction:

- Painted steel cabinet and aluminum drawers

#### Features:

- Dual column medication cart with 72 multi-dose bin capacity, single width narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Ball-bearing, full extension drawer slides
- Painted steel pull-out side shelf
- Plastic replaceable top with integrated push handles
- Corner bumpers, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL72MDSP) Included are: #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680435 aluminum waste container and #684810 locking sharps container



Capacity of 12 bins per drawer.  
(Bins not included)

#### Internal Drawer Dimensions:

##### Left Column Drawers:

- 6 ea – 4" drawers with internal dividers, capacity 12 bins per drawer (Bins are not included)
- 1 ea – 9.75" multipurpose storage drawer with adjustable dividers

##### Right Column Drawers:

- 2 ea – 4" multipurpose storage drawer
- 1 ea – 8" locking narcotics drawer
- 2 ea – 8" multipurpose storage drawer

#### Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>

#### Assembled Cart Dimensions:

- 45.75" H x 43" W x 23.5" D  
(116.2 cm x 109.2 cm x 59.7 cm)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 210 lbs (95.5 kg)
- Approx shipping weight: 245 lbs (111.13 kg)
- Approx Shipping Dimensions: 57" H x 47" W x 27" D