

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

V-Series, Five Drawer Procedure Cart with Key Lock #V30-5K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.



Sixteen available color choices: *Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Specifications

Construction:

Steel construction

Features:

- Key lock with two keys
- Aluminum side rail / push handle
- 30" of vertical drawer space
- Stabilizing bumper frame with extended front casters and
 (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

Internal Drawer Dimensions:

• 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D (101.6cm H x 74.9 W x 55.9 D)
- Cart weight: 115 lbs (52.2 kg)
- Approx shipping weight: 138 pounds (62.6 kg)

Drawer Configuration:

- (2) 3" (7.6 cm)
- (1) 6" (15.2 cm)
- (2) 9" (22.9 cm)

Integrated Pull Handle on all drawers