



650 Ford Street  
Colorado Springs, CO 80915  
Toll: (800)433-4064 \* Ph: (719)637-0300  
Fax: (719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## Value Line Med Bin Cart with 12 – 5 inch Bins, Basic Electronic Pushbutton Lock with Key Lock Override #VLT12EBIN5

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Basic Electronic Pushbutton Lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- Internal locking narcotics box in bottom 9" drawer
- Pull out side shelf
- Includes 12 Harloff 5" medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Certified "Latex Free"
- Five year warranty
- Compatible with most accessories
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish



Color Shown is Sand.



#### Thirteen available color choices:

(may be viewed at [www.harloff.com](http://www.harloff.com))

- Sand, Gloss Beige, Burgundy, Light Blue, Yellow, Navy, Light Gray, Mauve, Red, Teal, Hammer Tone Blue, Hammer Tone Gray and Hammer Tone Green

#### Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 3 - ea 3" bin drawers including 4 - 5" bins each
- 1 - ea 9" drawer with 2 adjustable dividers
- 1 - ea 9" drawer with 1 divider and locking narcotics box

#### Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

#### Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D  
(121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 145 pounds (65.9 kg)
- Approx. shipping weight: 185 pounds (81.8 kg)
- Approx. shipping dimensions: 53" H x 33" W x 26" D