

Value Line Medication Cart with Key Lock (MEDICINE ON TIME® compatible) #VLT150MOT

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system with 4 keys
- Capacity of up to 150 Medicine On Time® cards
- Internal locking narcotics box made for M.O.T. packaging
- Pull out side shelf
- Includes 20 Harloff 5" medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawers
- 1 - ea 9.75" drawer with locking narcotics compartment made for M.O.T. packaging.
- 2 - ea 9.75" drawers with 1 adjustable divider

Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D
(121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 140 pounds (63.7 kg)
- Approx. shipping weight: 180 pounds (81.8kg)

Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>