

Value Line Med Bin Cart with 18 – 3.5 inch Bins, Key Locking #VLT18BIN3

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system with 4 keys
- Internal narcotics box in top 9" drawer
- Pull out side shelf
- Includes 18 Harloff 3.5" medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Antimicrobial plastic top, removable
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 3 - ea 3" bin drawers including 6 -3.5" bins each
- 1 - ea 9" drawer with 1 divider and locking narcotics box
- 1 - ea 9" drawer with 2 adjustable dividers

Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D
(121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 155 pounds (70.5 kg)
- Approx. shipping weight: 195 pounds (86.6 kg)
- Approx. shipping dims: 57"H x 33" W x 26" D
(144.8cm H x 83.8cm W x 66cm D)

Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options: <https://www.harloff.com/finish-options/>

