



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Value Line Med Bin Cart with 20 – 5 inch Bins, Basic Electronic Pushbutton Lock with Key Lock Override
#VLT20EBIN5

Specifications



Color Shown is Beige



Construction:

- Painted steel cabinet and drawers

Features:

- Basic Electronic Pushbutton Lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- Internal narcotics box in bottom 9” drawer
- Pull out side shelf
- Includes 20 Harloff 5” medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5” ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five year warranty
- Compatible with most accessories
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish

Drawer Configuration:

- 2 - ea 3” multipurpose storage drawers
- 5 - ea 3” bin drawers including 4 - 5” bins each
- 1 - ea 9” drawer with 1 divider and locking narcotics box

Internal Drawer Dimensions:

- 23” W x 17” D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75” H x 29.5” W x 22.5” D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25” (109.9 cm)
- Cart weight: 172 pounds (78.2 kg)
- Approx. shipping weight: 207 pounds (94.1kg)

Thirteen available color choices:

(may be viewed at www.harloff.com)

- Sand, Gloss Beige, Burgundy, Light Blue, Yellow, Navy, Light Gray, Mauve, Red, Teal, Hammer Tone Blue, Hammer Tone Gray and Hammer Tone Green