



# MEDICATION CARTS

**VLT360EPC**

## VALUE PUNCH CARD MED CART

450 Punch Card Capacity with an E-Lock

### Features:

- Single column medication cart for punch card and liquid storage with a maximum capacity of 450 punch cards
- Includes one in-drawer single width locking narcotics box in top 9" drawer; comes with an Illinois lock with four keys
- Basic Electronic Pushbutton Lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- Painted steel cabinet with aluminum drawers
- Durable scratch resistant heat sealed powder coat finish
- Pull-out shelf
- Taller drawer height for increased capacity
- Ball-bearing, full extension drawer slides
- Painted steel top rail/push handle
- Removable plastic top
- Protective gray corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Adjustable punch card row dividers
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Certified "Latex Free"
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Color shown is Hammertone Green. More than 20 additional colors are available.



### Value 450 Punch Card Capacity Medication Cart with Electronic Pushbutton Lock on Cabinet

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

#### Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 1 - ea 9.75" drawer with 1 divider and locking narcotics box
- 2 - ea 9.75" drawers with two adjustable dividers

#### Assembled Cart Dimensions:

- Writing surface: 43.25" (109.9 cm)
- 46.75" H x 29.5" W x 22.5" D (121.3 cm H x 75 cm W x 57.2 cm D)
- Cart weight: 140 pounds (63.7 kg)
- Approx shipping weight: 175 pounds (79.4 kg)
- Approx shipping dims: 56"H x 33" W x 26" D (144.8cm H x 83.8cm W x 66cm D)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Engineered & Manufactured in the U.S.**



View similar products online at

<https://www.harloff.com/product-category/punch-card-carts-cabinets/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)**  
**[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**