

Toll: (800) 433-4064 Ph: (719) 637-0300 dicalsales@harloff.com

medicalsales@harloff.com www.harloff.com

Value Line Punch Card Medication Cart with Key Lock #VLT360PC

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

https://www.harloff.com/finish-options/

Construction:

· Painted steel cabinet and drawers

Features:

- · Tall cabinet with key locking system with 4 keys
- · Capacity of up to 360 punch card carts
- Internal locking narcotics box in top 9" drawer
- · Painted steel pull-out side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel, one with brake, one directional
- Painted steel top rail/push handle
- Antimicrobial plastic top, removable
- Full extension ball-bearing drawer slides
- · Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- · Certified "Latex Free"
- Twelve year manufacturer's warranty

Internal Drawer Dimensions:

• 23" W x 17" D (58.5 cm W x 43.2 cm D)

Drawer Configuration:

- 1 ea 3" multipurpose storage drawer
- 1 ea 9.75" drawer with 1 divider and locking narcotics box
- 2 ea 9.75" drawers with two adjustable dividers

Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D
 (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 140 pounds (63.7 kg)
- Approx. shipping weight: 170 pounds (81 kg)
- Approx. shipping dims: 56"H x 33" W x 26" D (144.8cm H x 83.8cm W x 66cm D)