



Wood Vinyl

WV600PC-DG

## PUNCH CARD MEDICATION CART

600 Capacity Wood Vinyl Aluminum, Key Lock

### Features:

- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- Double column medication cart with 600 punch card capacity
- Key locking cabinet with locking narcotics drawer, both with 4 keys
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- Pull-out shelf
- 4" ball-bearing premium rubber casters under skirting
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture
- Included skirt option to cover wheels, #WV-4WSKIRT-DG

### Colors:

Shown in Driftwood Gray color. Harloff offers three standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



Driftwood Gray (DG)



Engineered & Manufactured in the U.S.



#### Drawer Configuration:

- LEFT COLUMN: One 4" multipurpose storage drawer and three 9.75" punch card drawers
- RIGHT COLUMN: One 4" multipurpose storage drawer, one 9.75" locking narcotics drawer and two 9.75" multipurpose storage drawers

#### Overall Cart Dimensions:

- 44" H x 44.75" W x 25.25" D (111.8 cm H x 113.6 cm W x 64.8 cm D)
- Writing surface: 43.5" (110.5 cm)
- Cart weight: 167 pounds (75.9 kg)
- Approx. shipping weight: 230 pounds (104.5 kg)
- Approx. shipping dims: 51" H x 52" W x 32" D (129.54 cm H x 132.08 cm W x 81.28 cm D)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.

View similar products online at

<https://www.harloff.com/product-category/medication-furniture-carts-cabinets/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)

[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300

